

---

## **Phase Visual Meter and algorithmic estimation (PVMae) for LISA**

---

### **Lab Manual: Setup and Linux PC**

*Romi Gilat, Joshua Marmor, Carlos Arias*  
*Laboratory of Space Systems and Optomechanics*

## **PHYSICAL SETUP**

*Turn on the following instruments on the laboratory table:*

1. Laser - *small gray button on the side, set to 190 mA*
2. Oscilloscope - *power button on bottom left*
3. Moku - *logo is the power button*
4. Power Supply - *set to 15 V = 0.01 A*
5. Photodetector - *small black switch on the side*

*(The CCD Camera should be initially connected to the computer through the USB)*

### **Opening vivado**

1. In terminal, type  
source /home/Xilinx/Vivado/2022.1/settings64.sh
2. Type “vivado”

### **OPEN GUI**

#### **(1): Connect to the Red Pitaya**

On the linux terminal, type the following command:

```
ssh root@<ip address of RP>
```

In our case, the ip address of the RedPitaya in use (**rp-f084ec**) is: **169.254.2.239**

When asked for a password, use **‘root’**

Now you are connected to the terminal on the linux OS inside the RedPitaya.

*NOTE: If you are unable to connect to the RedPitaya, you may need to reboot the computer using **dzdo reboot**. Sometimes, when the system has been left alone for a while, the connection between the two times out.*

#### **(2): Navigate to the folder with the bitstream file**

Now that you are on the RedPitaya’s terminal, run the following command:

```
cd PVMAE/
```

Inside this folder, there are the executable files.

Any .bit —> compiles FPGA

Any .bpf —> communications between the RP and the Linux PC

#### **(3): Compile the FPGA bitstream**

Still in the RedPitaya’s terminal, type the following command:

```
cat pvmae_Test_100us.bit > /dev/xdevcfg  
{compile} {      what      } {    where    }
```

#### **(4): Initialize communication between the RedPitaya and the Linux PC**

In the RedPitaya terminal, run the following command:

```
./bpf_server.out
```

This will activate the communication protocol between the PC and the RedPitaya

#### **(5): Open a new terminal window in the folder with the Linux Code**

Navigate to the following folder address:

```
LASSO-DAQ/PVMae/Alberto
```

It makes things easier if you just access the folder, right click, and open a terminal from the menu.

#### **(6): Ensure all codes are compiled correctly**

In the terminal window opened in step (5), run this command:

```
make
```

If any changes have been made, this command will compile those new changes to make sure all files are up to date.

#### **(7): Connect the camera and FPGA and open the GUI**

In the same window, run the following command:

```
dzdo ./PVMae
```

Now the PC is connected to the RedPitaya and the CCD camera. The PC is the *master* and both the RP and the camera are *slaves*. This command will ask for a user password first, then the GUI will open.

### **GUI COMMANDS**

*To use these callbacks at least one window of the displays made by the program has to be active.*

→ 'q' : closes the camera, clears the memory, and closes the program.

- 's' : calls the function `save_data()` to save the current displayed information.
- '+' : changes the color palette for the phase front by increasing the contribution (range) of white color to scale and display the data.
- '-' : changes the color palette for the phase front by decreasing the contribution (range) of white color to scale and display the data.
- 'l' : increases the brightness on the phase front display.
- '\*' : decreases the brightness on the phase front display.
- 'b' : calls the function `back_gnd()` to take pictures of the background, averages it, and saves it in the matrix `background[i][j]`.
- 'u' : Increases the scale's upper limit for the phase front
- 'i' : Decreases the scale's upper limit for the phase front
- 'l' : Increases the scale's lower limit for the phase front
- 'k' : Decreases the scale's lower limit for the phase front
- 'h' : Increases scale's range for the phase front
- 'j' : Decreases scale's range for the phase front

## **CONNECTING TO THE MOKU**

### **(1): Power on the Moku**


Press the logo in the middle to power on

### **(2): Connect on the iPad:**

Password: 111111

Wifi: Moku-000927

### **(3): Open Oscilloscope App**



### **LINUX PC INFO**

ID: **aero-5w8bmh3**

### **REDPITAYA INFORMATION**

*Primary RP*

ID: **rp-f084ec**

IP: **169.254.2.239**

*Secondary RP*

ID: **rp-f0949f**

IP: **169.254.8.101**

### **FLYCAPTURE**

Software: *FlyCap2*

Camera Serial Number: *Chameleon CMLN-13S2M*

**NOTE:** *To open the software, use run the following command:*

```
dzdo flycap
```

Running the command without dzdo will not detect the camera.